

Boule Cal

CAL

IVD



Open vial stability 5 days ¹

Calibrator

Lot	
21806-34	2018-09-07

Medonic

Parameter ²	RBC 10 ¹² /L	MCV fL	Plt 10 ⁹ /L	MPV fL	WBC 10 ⁹ /L	Hgb g/dL
CA620/530	3,98	76,0	260	10,6	9,0	11,3
Medonic M-series / M32	3,98	78,9	240	10,6	8,7	11,4
Range ± ³	0,10	2,0	15	1,5	0,3	0,2

Swelab

Parameter ²	RBC 10 ¹² /L	MCV fL	Plt 10 ⁹ /L	MPV fL	WBC 10 ⁹ /L	Hgb g/dL
AC9xxeo+	3,98	79,8	233	11,5	7,9	11,3
Swelab Alfa / Alfa Plus	3,98	79,1	240	11,2	8,7	11,4
Range ± ³	0,10	2,0	15	1,5	0,3	0,2

	1	2	3
bg	Трайност на отворен флакон	5 дни	Обхват
hr	Stabilnost otvorene bočice	5 dani	Parametar
cs	Stabilita otevřené ampulky	5 dny	Parametr
da	Åben rør stabilitet	5 dage	Parameter
et	Stabiilsus avatud vialis	5 päevade arv	Parameeter
fr	Stabilité en flacon ouvert	5 jour	Paramètre
de	Stabilität geöffneter Flaschen	5 tage	Parameter
el	Διατήρηση δείγματος μετά την αποσφράγιση	5 ημέρα	Παράμετρος
it	Stabilità della fiala aperta	5 giorni	Parametro
lv	Atvērt flakonu stabilitāti	5 dienas	Parametrs
lt	Stabilumas atidarius buteliuka	5 d.	Parametras
no	Åpen rørglass stabilitet	5 dager	Parameter
pl	Trwałość po otwarciu fiołki	5 Liczba dni	Kryterium
pt	Estabilidade após abertura do frasco	5 dias	Parâmetro
ro	Termenul de valabilitate al fiolei desfăcute	5 zile	Parametru
ru	Стабильность вскрытого флакона	5 дней	Параметр
sr	Proizvod stabilan nakon otvaranja pakovanja	5 dana	Parametar
sk	Stabilita otvorenej injekčnej liekovky	5 dní	Parameter
sl	Stabilnost odprte stekleničke	5 dni	Parameter
es	Estabilidad de la cápsula abierta	5 dias	Parámetro
sv	Hållbarhet för öppnad flaska	5 dagar	Parameter
tr	Açık şişe stabilitesi	5 gün	Parametre



<http://www.medonic.se/>
www.swelab.com/


Boule Medical AB
Domnarvsgatan 4
SE-163 53 Spånga, Sweden

25123-3

1504055

For CA620 with barcode reader

1. Go to Menu 7.4.1

<p>C.A620 Main Menu Program: V3c75</p> <ol style="list-style-type: none"> 1. Sample memory... 2. Set auto-print-mode=0 - 3. Set print-key function=2 - 4. Set print format= 8 - 5. Setup menu... 6. Service menu... 7. QC menu... 8. System flow menu... 	<p>7. QC menu</p> <ol style="list-style-type: none"> 7.1 Control memory 7.2 Control L-J plots 7.3 X-B L-J plots 7.4 QC setup menu... 7.5 Set SEQ= 1 	<p>7.4.1 Input control definition</p> <p>Use the barcode reader to input the list of barcodes from the assay sheet for the control you want to define. See User's Manual for further info.</p> <p>0 barcodes read. Press <Menu> to skip.</p>
---	---	--

2. Scan in the calibrator/control in sequence top to bottom

Mfg: CA Boule reag.

Product: Cal

Lot: 21806-34

ID: 2180634+

Exp.: 7-Sep-2018



For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 21806-34
Prod.: Cal
ID: 2180634+
Exp.: 7-Sep-2018



For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa

Lot: 21806-34

Prod.: Cal

ID: 2180634+

Exp.: 7-Sep-2018

