

CONTROL

Boule 3PD Hematology Control

IVD

РУ ФСЗ 2008/02770 от 21 Июля 2016 г.

Стабильность открытого флакона: 14 дней

4,5 мл

| | | |
|------------------|-----------------|-------------------|
| Tri-Level | 31809-5K | 2019.02.06 |
| Low | 31809-51 | 2019.02.08 |
| Normal | 31809-52 | 2019.02.06 |
| High | 31809-53 | 2019.02.09 |

| | | CONTROL | Low |
|--------------------------|---------------------|-------------------|----------------------|
| | | LOT | 31809-51 |
| | | Mean ² | ± Range ³ |
| SWELAB AC 9xx eo+ | | | |
| RBC | 10 ¹² /L | 2,27 | 0,12 |
| MCV | fL | 69,2 | 7,0 |
| RDW | % | 14,3 | 4,5 |
| Hct | % | 15,6 | 3,0 |
| Plt | 10 ⁹ /L | 89 | 15 |
| MPV | fL | 11,3 | 1,5 |
| WBC | 10 ⁹ /L | 3,2 | 0,3 |
| Hgb | g/dL | 5,6 | 0,3 |
| MCH | pg | 24,7 | 2,0 |
| MCHC | g/dL | 36,0 | 4,0 |
| Lymph# | 10 ⁹ /L | 1,4 | 0,4 |
| Mid# | 10 ⁹ /L | 0,3 | 0,3 |
| Gran# | 10 ⁹ /L | 1,4 | 0,8 |
| Lymph | % | 41,5 | 7,0 |
| Mid | % | 13,9 | 5,0 |
| Gran | % | 44,6 | 7,0 |

| | | CONTROL | Normal |
|--|--|-------------------|----------------------|
| | | LOT | 31809-52 |
| | | Mean ² | ± Range ³ |
| | | 4,01 | 0,18 |
| | | 79,5 | 7,0 |
| | | 11,8 | 4,5 |
| | | 32,1 | 3,5 |
| | | 257 | 30 |
| | | 12,0 | 1,5 |
| | | 8,2 | 0,6 |
| | | 11,3 | 0,4 |
| | | 28,2 | 2,0 |
| | | 35,1 | 4,0 |
| | | 4,0 | 0,6 |
| | | 0,6 | 0,5 |
| | | 3,5 | 1,0 |
| | | 44,9 | 5,0 |
| | | 7,0 | 5,0 |
| | | 48,1 | 5,0 |

| | | CONTROL | High |
|--|--|-------------------|----------------------|
| | | LOT | 31809-53 |
| | | Mean ² | ± Range ³ |
| | | 4,69 | 0,22 |
| | | 89,1 | 7,0 |
| | | 12,6 | 4,5 |
| | | 41,8 | 4,0 |
| | | 541 | 60 |
| | | 12,9 | 1,5 |
| | | 19,4 | 1,8 |
| | | 15,0 | 0,5 |
| | | 31,9 | 2,0 |
| | | 35,9 | 4,0 |
| | | 12,1 | 1,5 |
| | | 1,3 | 1,0 |
| | | 6,0 | 2,5 |
| | | 62,0 | 5,0 |
| | | 6,7 | 5,0 |
| | | 30,8 | 5,0 |

