

CONTROL

Boule 3PD Hematology Control

IVD

РУ ФСЗ 2008/02770 от 21 Июля 2016 г.

Стабильность открытого флакона: 14 дней

Tri-Level	31809-8K	2019.02.06
Low	31809-51	2019.02.08
Normal	31809-54	2019.02.06
High	31809-53	2019.02.09

		CONTROL	Low
		LOT	31809-51
SWELAB AC 9xx eo+ Medor		Mean ²	± Range ³
RBC	10 ¹² /L	2,27	0,12
MCV	fL	69,2	7,0
RDW	%	14,3	4,5
Hct	%	15,6	3,0
Plt	10 ⁹ /L	89	15
MPV	fL	11,3	1,5
WBC	10 ⁹ /L	3,2	0,3
Hgb	g/dL	5,6	0,3
MCH	pg	24,7	2,0
MCHC	g/dL	36,0	4,0
Lymph#	10 ⁹ /L	1,4	0,4
Mid#	10 ⁹ /L	0,3	0,3
Gran#	10 ⁹ /L	1,4	0,8
Lymph	%	41,5	7,0
Mid	%	13,9	5,0
Gran	%	44,6	7,0

		CONTROL	Normal
		LOT	31809-54
		Mean ²	± Range ³
		3,97	0,18
		79,5	7,0
		12,8	4,5
		31,8	3,5
		259	30
		12,0	1,5
		8,0	0,6
		11,1	0,4
		27,9	2,0
		34,8	4,0
		3,9	0,6
		0,6	0,5
		3,4	1,0
		44,9	5,0
		7,0	5,0
		48,1	5,0

		CONTROL	High
		LOT	31809-53
		Mean ²	± Range ³
		4,69	0,22
		89,1	7,0
		12,6	4,5
		41,8	4,0
		541	60
		12,9	1,5
		19,4	1,8
		15,0	0,5
		31,9	2,0
		35,9	4,0
		12,1	1,5
		1,3	1,0
		6,0	2,5
		62,0	5,0
		6,7	5,0
		30,8	5,0

