

Boule 3PD Hematology Control

CONTROL

IVD

РУ ФСЗ 2008/02770 от 21 Июля 2016 г.

4.5 mL

Стабильность открытого флакона: 14 дней

Tri-Level		31803-5K	2018.08.07
Low		31803-51	2018.08.09
Normal		31803-52	2018.08.07
High		31803-53	2018.08.07

LOW CONTROL		Lot 31803-51	Sysmex KX-21, KX-21N, K-800, K-1000, XP-100, XP-300, почH-100i	Mean	± Range
Parameters					
WBC	x10 ⁹ /L			3,2	0,3
Lymph #	x10 ⁹ /L			1,5	0,5
Mid # (Mxd*)	x10 ⁹ /L			0,3	0,3
Gran # (Neut*)	x10 ⁹ /L			1,5	0,8
Lymph	%			47,2	5,0
Mid (Mxd*)	%			5,5	4,0
Gran (Neut*)	%			47,3	5,0
RBC	x10 ¹² /L			2,18	0,12
Hgb	g/dL			4,9	0,3
Hct	%			14,7	3,0
MCV	fL			67,5	7,0
MCH	pg			22,5	2,0
MCHC	g/dL			33,3	4,0
RDW-CV	%			12,2	4,5
RDW-SD	fL			33,7	5,7
PLT	x10 ⁹ /L			73	15
MPV	fL			9,0	1,5
NORMAL CONTROL		Lot 31803-52	Sysmex KX-21, KX-21N, K-800, K-1000, XP-100, XP-300, почH-100i	Mean	± Range
Parameters					
WBC	x10 ⁹ /L			8,8	0,6
Lymph #	x10 ⁹ /L			4,3	0,6
Mid # (Mxd*)	x10 ⁹ /L			0,8	0,5
Gran # (Neut*)	x10 ⁹ /L			3,7	1,0
Lymph	%			48,8	5,0
Mid (Mxd*)	%			9,6	4,0
Gran (Neut*)	%			41,6	5,0
RBC	x10 ¹² /L			4,04	0,18
Hgb	g/dL			11,7	0,4
Hct	%			31,9	3,5
MCV	fL			79,0	7,0
MCH	pg			29,0	2,0
MCHC	g/dL			36,7	4,0
RDW-CV	%			10,4	4,5
RDW-SD	fL			34,5	5,2
PLT	x10 ⁹ /L			240	30
MPV	fL			9,1	1,5
HIGH CONTROL		Lot 31803-53	Sysmex KX-21, KX-21N, K-800, K-1000, XP-100, XP-300, почH-100i	Mean	± Range
Parameters					
WBC	x10 ⁹ /L			20,9	1,8
Lymph #	x10 ⁹ /L			13,0	1,5
Mid # (Mxd*)	x10 ⁹ /L			1,8	1,0
Gran # (Neut*)	x10 ⁹ /L			6,1	2,5
Lymph	%			62,0	5,0
Mid (Mxd*)	%			8,8	4,0
Gran (Neut*)	%			29,2	5,0
RBC	x10 ¹² /L			4,69	0,22
Hgb	g/dL			15,3	0,5
Hct	%			40,6	4,0
MCV	fL			86,5	7,0
MCH	pg			32,6	2,0
MCHC	g/dL			37,7	4,0
RDW-CV	%			11,2	4,5
RDW-SD	fL			37,7	4,8
PLT	x10 ⁹ /L			525	60
MPV	fL			9,4	1,5